


# Deep, Unbiased and Quantitative mass spectrometry-based plasma proteome analyses of adaptive response to COVID-19 vaccine



Ting Huang<sup>\*1</sup>, Alex Rosa Campos<sup>2</sup>, Ramón Díaz Peña<sup>4</sup>, Svetlana Maurya<sup>4</sup>, Khatereh Motamedchaboki<sup>1</sup>, Jian Wang<sup>1</sup>, Laura R. S. Oliveira<sup>3</sup>, Camila I. Alves<sup>3</sup>, Rafael L. Rosa<sup>3</sup>, Yohana P. C. Alves<sup>3</sup>, Lucélia Santi<sup>3</sup>, and Walter O. Beys-da-Silva<sup>3</sup>


<sup>1</sup>Seer, Inc., Redwood City, CA 94065, USA; <sup>2</sup>Plexium, San Diego, CA, USA; <sup>3</sup>Federal University of Rio Grande do Sul, Porto Alegre, Brazil; <sup>4</sup>Sanford Burnham Prebys, La Jolla, CA, USA.

## Proteograph™ workflow coupled with TMTpro 18Plex provides deep access to plasma proteome

## Proteograph workflow provide valuable insights on personalized human responses to Covid-19 vaccination

### Introduction


COVID-19 vaccines have been extensively used to immunize a large worldwide population during COVID-19 pandemic. However, individuals respond to the vaccine differently, leading to distinct proteome changes and vaccination efficiency. At the same time, our understanding of the underlying molecular mechanisms behind the adaptive response to COVID-19 vaccine is still limited. Here, we used the Proteograph workflow<sup>1</sup> coupled with TMT labeling to quantify the proteome changes in a cohort of 12 volunteers immunized with two doses of COVID-19 vaccine. Our data confirmed that Proteograph workflow enables identification of more than 3,000 proteins from human plasma samples in this study. Since 7 volunteers still got COVID-19 after two-dose vaccine, our analysis gained further insights into the personalized response to vaccination.



### Methods: Proteograph Workflow and TMT Labeling

Plasma samples were collected before and after first and second vaccination from each volunteer. Samples were processed by Proteograph™ Product Suite (Seer Inc.) using five distinct nanoparticles (NPs). Tryptic peptides were then pooled into one single sample for TMT labeling. A total of 48 samples were allocated into three 18-plex TMT mixtures, followed by peptide fractionation by high pH RP and LC-MS/MS analysis, comprised of an EASY nanoLC system coupled to an Orbitrap Fusion Lumos MS equipped with FAIMS Pro Interface (all from Thermo Fisher Scientific). Peptides were separated using a 25-cm analytical C18 Aurora column (IonOpticks) into 24 fractions. A two-hour gradient with three FAIMS Pro compensation voltages was used in this study. The raw spectra data were processed by SpectraMine (Biognosys) to generate peptide identification and quantification data. We then used R/Bioconductor package MSstatsTMT for protein quantification, normalization, and differential analysis.


### A Proteograph Workflow with TMTpro 18plex Labeling



### B Experimental Design

| Plex | 126        | 127N        | 127C       | 128N        | 128C       | 129N        | 129C       | 130N        | 130C        | 131N         | 131C        | 132N         | 132C        | 133N         | 133C        | 134N         | 134C  | 135N  |
|------|------------|-------------|------------|-------------|------------|-------------|------------|-------------|-------------|--------------|-------------|--------------|-------------|--------------|-------------|--------------|-------|-------|
| set1 | S1 preVax1 | S1 postVax1 | S1 preVax2 | S1 postVax2 | S2 preVax1 | S2 postVax1 | S2 preVax2 | S2 postVax2 | S3 preVax1  | S3 postVax1  | S3 preVax2  | S3 postVax2  | S4 preVax1  | S4 postVax1  | S4 preVax2  | S4 postVax2  | pool1 | pool2 |
| set2 | S6 preVax1 | S6 postVax1 | S6 preVax2 | S6 postVax2 | S8 preVax1 | S8 postVax1 | S8 preVax2 | S8 postVax2 | S9 preVax1  | S9 postVax1  | S9 preVax2  | S9 postVax2  | S10 preVax1 | S10 postVax1 | S10 preVax2 | S10 postVax2 | pool1 | pool2 |
| set3 | S5 preVax1 | S5 postVax1 | S5 preVax2 | S5 postVax2 | S7 preVax1 | S7 postVax1 | S7 preVax2 | S7 postVax2 | S11 preVax1 | S11 postVax1 | S11 preVax2 | S11 postVax2 | S12 preVax1 | S12 postVax1 | S12 preVax2 | S12 postVax2 | pool1 | pool2 |

### C Waterfall Plot (HPPP) - Number of Quantified Protein Groups: 3094



### Figure 1. Proteograph workflow with TMT labeling: An overview of protein quantification performance.

A) For each biological sample, five peptide mixtures prepared by five distinct nanoparticles of Proteograph™ Assay were pooled into one single sample. Then, 16 pooled samples from 4 volunteers were labeled with 18plex TMT and combined to produce a single TMT mixture. B) 48 samples from 12 volunteers were allocated into three 18plex TMT mixtures. Each mixture includes two masterpool samples for normalization which was labeled with extra 2 labels in TMTpro 18plex labeling kit. C) Compared with the existing plasma proteomics studies on COVID-19 vaccination<sup>2,3</sup>, the combination of Proteograph workflow with TMTpro 18Plex has provided much deeper access to human plasma proteome, i.e., 23,372 peptides and 3,094 proteins. Mapping these 3,094 proteins to the Human Plasma Proteome Project (HPPP) protein database<sup>4</sup> shows detection across the entire concentration range of the database, 20% of which fell in the low abundance range.




Figure 2. Clustering and differential protein analysis across the cohort.

A) Hierarchical clustering analysis showed that samples from the same volunteer were clustered together and demonstrated a high inter-individual variation. B) Differential analysis using samples from 12 volunteers before and after two doses of COVID-19 vaccine. None of the proteins were differentially expressed after the first dose but 69 proteins became differential after the second dose.




Figure 3. Individual-specific differential proteomic analysis.

A) Proteins from each volunteer were clustered into significant discrete groups using mFuzz package to illustrate various patterns of protein expression changes after vaccination. Most volunteers shared a common N (up-regulation) pattern, where protein abundances were up-regulated by the first dose, decreased after the first dose, and boosted by the second dose. B) Pathway enrichment analyses on individual volunteers showed that individuals who enriched most GO-BP and immune pathways didn't get COVID-19 within 10 months after two-dose vaccination. C) Clustering on immune pathways-related proteins, individuals who didn't get Covid-19 were grouped together.

### Conclusion

- Proteograph workflow coupled with TMT labeling quantifies >3,000 protein groups from a cohort of COVID-19 vaccination study with 48 plasma samples
- Covid-19 vaccine related plasma proteomics data show a high between-subject variance
- Individuals with enriched gene oncology and immune pathways didn't get covid within 10 months after vaccination

### References

- Blume et al. Nat. Comm. (2020)
- Wang et al. Front. Immunol. (2022)
- Wang et al. medRxiv (2022)
- Schwenk, et al. Journal of Proteome Research. (2017)



Publications

Copyright Seer, Inc 2023

Seer, Inc., Redwood City, CA 94065, USA | \*tuang@seer.bio

